

DOSAGE CHART · CANITIDE REPAIR

BPC-157 for dogs: dose by weight

BPC-157 10 mg + TB-500 10 mg. 5 days on / 2 days off · 8-week cycle, then 4 weeks off.

WEIGHT BAND	PER DOSE	IN MICROGRAMS	PER WEEK	A CARTRIDGE LASTS
Small under 20 lb	6 units	198 mcg	30 units	~10 weeks
Medium 20–50 lb	9 units	297 mcg	45 units	~6 weeks
Large 50–90 lb	12 units	396 mcg	60 units	~5 weeks
Giant 90 lb and up	15 units	495 mcg	75 units	~4 weeks

Reading the chart. One pen unit is 0.01 mL and carries about 33 mcg of BPC-157. A 10 mg vial of BPC-157 reconstituted in 3 mL of bacteriostatic water is 3.33 mg/mL, the same concentration as a REPAIR cartridge, so the unit numbers above read the same on a U-100 insulin syringe.

Missed a dose? Skip it and continue on schedule. Never double up.

Your vet first. Bring this chart to your veterinarian before starting, especially if your dog takes an NSAID, a steroid or insulin.

Dog's name

Weight and band

Start date