



THE CANITIDE DOSING & PEN GUIDE

Your dog's numbers. Ten seconds a day.

Weight-band dosing for all six Canitide formulas, how to use the pen, how to store cartridges, what to watch for, and when to call your vet. Keep it with the pen.

Every Canitide formula is dosed in pen units per weight band. One pen unit is 0.01 mL. Cartridges hold 300 usable units after the 2-unit prime; they are filled to 3.15 mL so priming never costs a dose. Find your dog's band, then the formula. Bring this card to your veterinarian.

The four bands

BAND	WEIGHT
Small	under 20 lb
Medium	20–50 lb
Large	50–90 lb
Giant	90 lb and up

Weigh your dog on the same scale each month. A band change changes the dose.

Reading the card

Per dose is what you dial. **Per week** adds the doses in a week. **Cartridge lasts** is how long one 300-unit cartridge covers at that band; the Club cadence is set so the next one arrives before this one runs out.

Missed a dose? Skip it and continue. Never double up.

Dosing by formula

● **Canitide REPAIR** · BPC-157 10 mg + TB-500 10 mg · 5 days on / 2 days off · 8-week cycle, then 4 weeks off

BAND	PER DOSE	PER WEEK	CARTRIDGE LASTS	CLUB CADENCE
Small under 20 lb	6 units	30 units	~10 weeks	every 8 wk
Medium 20–50 lb	9 units	45 units	~6 weeks	every 6 wk
Large 50–90 lb	12 units	60 units	~5 weeks	every 4 wk
Giant 90 lb and up	15 units	75 units	~4 weeks	every 4 wk

● **Canitide COAT** · GHK-Cu 50 mg · Daily · 8-week cycle, then maintenance

BAND	PER DOSE	PER WEEK	CARTRIDGE LASTS	CLUB CADENCE
Small under 20 lb	4 units	28 units	~10 weeks	every 8 wk
Medium 20–50 lb	6 units	42 units	~7 weeks	every 6 wk
Large 50–90 lb	8 units	56 units	~5 weeks	every 4 wk
Giant 90 lb and up	10 units	70 units	~4 weeks	every 4 wk

● **Canitide VITAL** · MOTS-c 10 mg · 3× per week · 8-week cycle

BAND	PER DOSE	PER WEEK	CARTRIDGE LASTS	CLUB CADENCE
Small under 20 lb	10 units	30 units	~10 weeks	every 8 wk
Medium 20–50 lb	15 units	45 units	~6 weeks	every 6 wk
Large 50–90 lb	20 units	60 units	~5 weeks	every 4 wk
Giant 90 lb and up	25 units	75 units	~4 weeks	every 4 wk

● **Canitide LEAN** · AOD-9604 5 mg · Daily (morning, before food) · 12-week cycle

BAND	PER DOSE	PER WEEK	CARTRIDGE LASTS	CLUB CADENCE
Small under 20 lb	4 units	28 units	~10 weeks	every 8 wk
Medium 20–50 lb	6 units	42 units	~7 weeks	every 6 wk
Large 50–90 lb	8 units	56 units	~5 weeks	every 4 wk
Giant 90 lb and up	10 units	70 units	~4 weeks	every 4 wk

● **Canitide LONG** · Epithalon 10 mg · Daily for 10 days · repeat every 6 months

BAND	PER DOSE	PER WEEK	CARTRIDGE LASTS	CLUB CADENCE
Small under 20 lb	15 units	105 units	one 10-day cycle	every 26 wk
Medium 20–50 lb	20 units	140 units	one 10-day cycle	every 26 wk
Large 50–90 lb	25 units	175 units	one 10-day cycle	every 26 wk
Giant 90 lb and up	30 units	210 units	one 10-day cycle	every 26 wk

● **Canitide GUT** · BPC-157 (arginate) 500 mcg + KPV 250 mcg per capsule · 60 delayed-release capsules · 1 capsule daily for small dogs · 2 daily for medium dogs and up

BAND	PER DOSE	PER WEEK	BOTTLE LASTS	CLUB CADENCE
Small under 20 lb	1 capsule	7 capsules	~8 weeks	every 8 wk
Medium 20–50 lb	2 capsules	14 capsules	~4 weeks	every 4 wk
Large 50–90 lb	2 capsules	14 capsules	~4 weeks	every 4 wk
Giant 90 lb and up	2 capsules	14 capsules	~4 weeks	every 4 wk

Using the pen

1. **Load.** Unscrew the body, drop in a cartridge, screw closed.
2. **Attach.** Twist a fresh 32G pen needle onto the tip. Remove both caps.
3. **Prime.** Dial 2 units, point up, press until a drop appears. Once per new cartridge.
4. **Dial.** Turn the dial to your dog's units. Big numerals, a click per unit.
5. **Dose.** Lift loose skin at the scruff or flank, press the button, hold five seconds.
6. **Done.** Cap the needle, drop it in the sharps container, cartridge back in the fridge.

Where

Loose skin at the scruff (between the shoulder blades) or the flank. Rotate sites. Most dogs barely notice a 32G × 4 mm pen needle; doing it at breakfast helps.

Needles

One needle per dose. Never share a pen between dogs. Used needles go in the mail-back sharps container; when it reaches the fill line, seal it and use the prepaid label.

Storage

- Cartridges: keep refrigerated at 2–8 °C. In use, the loaded pen may stay at room temperature for up to 28 days.
- Capsules: room temperature, cap tight, out of the sun.
- Pen needles: single use. Never share a pen between dogs.

Running more than one formula

Protocols pair formulas on purpose, but start them two weeks apart so you can tell what is doing what. Use separate cartridges and doses; never mix cartridges. Keep a dated record of each dose and what you observe, ideally in Fieldnotes at canipep.com, so your vet can read it.

What to watch for

Stop and call your veterinarian if you see anything unusual after a dose: swelling beyond a small bump at the site, vomiting, lethargy that lasts past the first day or two, or any change that worries you. Then tell us at hello@canitides.com so we can log it.

Do not start any Canitide protocol for a dog with cancer or on chemotherapy without the oncologist's explicit sign-off. Not for pregnant or nursing dogs, or dogs under one year. If your dog takes NSAIDs, steroids or insulin, talk to your vet before starting.

Your dog's card

Dog	_____	Weight	_____ lb
Band	☐ Small ☐ Medium ☐ Large ☐ Giant	Vet	_____
Formula	_____	Units per dose	_____
Schedule	_____	Started	___ / ___ / ___

Cartridges are filled to 3.15 mL (315 units) so the 2-unit prime and hub loss never cost a dose; all dosing math uses 300 usable units per cartridge. Daily formulas are capped at 10 units per dose so a cartridge always covers a 4-week cadence.

Canitide products are sold for research purposes and are not FDA-approved drugs. They are not intended to diagnose, treat, cure or prevent any disease. No controlled canine clinical trials have been completed. Weight bands are calculated from the mg/kg ranges in the published animal research; they are a starting point for a conversation with a licensed veterinarian, not a prescription. Not for human use. Keep out of reach of children. © 2026 Canitides · canitides.com · hello@canitides.com